

# Panexin CD

*Serum Replacement  
with Defined Components*



## Ordering information

Panexin CD, Chemically defined serum replacement for adherent and non-adherent cells

**Cat. No. P04-93100 100 ml**  
**Cat. No. P04-930500 500 ml**

Phone: +49 (0) 8543 6016 -30  
Online Order at: [www.pan-biotech.de](http://www.pan-biotech.de)

## Key Product Features

- ✓ No more lot testing
- ✓ Chemically defined formulation
- ✓ High reproducibility and simplified downstream processing
- ✓ Ideal for the development of sensitive cell-based *in vitro* tests and coculture procedures

## More Efficiency




Figure 1: Efficiency and growth stimulation of Panexin CD compared to FBS (10 % in DMEM) for different cell types. Tests were performed at 37°C, 5% CO<sub>2</sub>. Cells were seeded with 10000/cm<sup>2</sup>



Visit our onlineshop at  
[www.pan-biotech.de](http://www.pan-biotech.de)

**PAN**  
BIOTECH™